

1. Part Number: 8111H-50S

2. Photograph:



3. Product Information:

Description:	SteriWare Powder Funnel 50mm - Sterile
Use:	Single use, disposable powder spatula & cover
Material of Construction:	High Density Polyethylene (HDPE) - virgin <ul style="list-style-type: none"> • FDA acceptable CFR 177.1520 • EU Regulations 10/2011 compliant • EC Regulations 1935/2004 compliant • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Max. Diameter:	68mm
ID at Top of Funnel:	50mm
Overall Height:	45mm
OD at Bottom of Spout:	18mm
ID at Bottom of Spout:	15mm
Nominal Weight:	5g
Moulding Environment:	Class 100,000 Cleanroom
Packing Environment	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25-45kGy)
Quantity per Box:	50
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture

ATEX:

This product has not been ATEX rated – Contact Sampling Systems for details

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable sampler. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 23-Jul-24