

1. Part Number: 8111P-350S

2. Photograph:



3. Product Information:

Description:	SteriWare® Powder Funnel - 350mm - Sterile
Use:	Single use, disposable funnel
Material of Construction:	Polypropylene (PP) - virgin <ul style="list-style-type: none"> • Conforms to FDA 21 CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Maximum Diameter:	343mm
Overall Height:	264mm
Outside Diameter of Tip of Spout:	52mm
Inside Diameter of Tip of Spout:	48mm
Length of Spout:	43mm
Nominal Weight:	348g
Moulding Environment:	General Factory
Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE bag)
Method of Sterilisation:	Gamma Irradiation (minimum dose 25 kGy)
Number of Funnels per Box:	8

Recommended Storage Conditions:

Dry and ambient temperature

Shelf Life:

2 years from the date of manufacture

ATEX:

This product has not been ATEX tested

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

change.control@sampling.com**5. Country of Origin:**

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable product.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 23-Jul-24