

**1. Part Number:** 8111P-350

**2. Photograph:**



### 3. Product Information:

Description:	SteriWare® Powder Funnel - 350mm
Use:	Single use, disposable funnel
Material of Construction:	Polypropylene (PP) - virgin <ul style="list-style-type: none"> <li>• Conforms to FDA 21 CFR 177.1520</li> <li>• Conforms to EU Regulations 10/2011</li> <li>• Conforms to EC Regulations 1935/2004</li> <li>• BSE/TSE free</li> </ul>
Method of Construction:	Single piece, injection moulded
Maximum Diameter:	343mm
Overall Height:	264mm
Outside Diameter of Tip of Spout:	52mm
Inside Diameter of Tip of Spout:	48mm
Length of Spout:	43mm
Nominal Weight:	348g
Moulding Environment:	General Factory
Primary Packing Environment:	Class 100,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE bag)
Method of Sterilisation:	N/A

Number of Funnels per Box: 8  
Recommended Storage Conditions: Dry and ambient temperature  
Shelf Life: 5 years from the date of manufacture  
ATEX: This product has not been ATEX tested

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

[change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable product.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 10-Sep-19