

**1. Part Numbers:** 8121P-4

**2. Photograph:**



**3. Product Information:**

Description	Narrow Mouth PP Bottle 4ml
Materials of Construction:	Polypropylene (PP) <ul style="list-style-type: none"> <li>• Conforms to FDA 21 CFR 177.1520</li> <li>• Conforms to USP Class VI</li> </ul>
Manufacturing Environment:	Controlled environment within factory
Nominal Volume:	4ml
Nominal Height:	41mm
Diameter:	15mm
Internal Neck Diameter:	8mm
Weight:	3g
Pack Size:	12 bottles per pack
Method of Sterilisation:	N/A (these bottles are non-sterile)
BSE/TSE Free:	Yes
Supplied with lid screwed on:	Yes
Recommended Storage Conditions:	Dry and ambient temperature

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 21-Jan-25