

1. Part Numbers: 8121P-4

2. Photograph:



3. Product Information:

Description	Narrow Mouth PP Bottle 4ml
Materials of Construction:	Polypropylene (PP) <ul style="list-style-type: none">• Conforms to FDA 21 CFR 177.1520• Conforms to USP Class VI
Manufacturing Environment:	Controlled environment within factory
Nominal Volume:	4ml
Nominal Height:	41mm
Diameter:	15mm
Internal Neck Diameter:	8mm
Weight:	3g
Pack Size:	12 bottles per pack
Method of Sterilisation:	N/A (these bottles are non-sterile)
BSE/TSE Free:	Yes
Supplied with lid screwed on:	Yes
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:

Name: Kuldeep Flora

Position: Quality Manager

Date: 23-Mar-22