

**1. Part Numbers:** 8126H-1000A

**2. Photograph:**



### 3. Product Information:

Description	Wide Mouth HDPE Bottle 1000ml Amber
Materials of Construction:	Bottle HDPE (conforms to FDA 177.1520 & USP Class VI) Cap PP (conforms to FDA 177.1520 & USP Class VI)
Manufacturing Environment:	Controlled environment within factory
Nominal Volume:	1000ml
Brimful Volume:	1075ml
Nominal Height:	201mm
Diameter:	92mm
Internal Neck Diameter:	52mm
Weight:	102g
Pack Size:	24 bottles
Method of Sterilisation:	N/A (these bottles are non-sterile)
BSE/TSE Free:	Yes
Supplied with lid screwed on:	Yes
Recommended Storage Conditions:	Dry and ambient temperature

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 2-Jul-21