

Specification Sheet

Specification Number: 8126H-30A
Revision: 05
Date Issued: 2-Jul-21

1. Part Numbers: 8126H-30A

2. Photograph:



3. Product Information:

Description Wide Mouth HDPE Bottle 30ml Amber

Materials of Construction: Bottle HDPE (conforms to FDA 177.1520 & USP Class VI)

Cap PP (conforms to FDA 177.1520 & USP Class VI)

Manufacturing Environment: Controlled environment within factory

Nominal Volume: 30ml
Brimful Volume: 34ml
Nominal Height: 63mm
Diameter: 37mm
Internal Neck Diameter: 21mm
Weight: 13g

Pack Size: 72 bottles

Method of Sterilisation: N/A (these bottles are non-sterile)

BSE/TSE Free: Yes
Supplied with lid screwed on: Yes

Recommended Storage Conditions: Dry and ambient temperature

4. Approval: Signed:

Name: Kuldeep Flora
Position: Quality Manager

Date: 2-Jul-21



4-5 Forge Mills Park, Coleshill. Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com