

**1. Part Numbers:** 8126H-60A

**2. Photograph:**



**3. Product Information:**

Description	Wide Mouth HDPE Bottle 60ml Amber
Materials of Construction:	Bottle HDPE (conforms to FDA 177.1520 & USP Class VI) Cap PP (conforms to FDA 177.1520 & USP Class VI)
Manufacturing Environment:	Controlled environment within factory
Nominal Volume:	60ml
Brimful Volume:	68ml
Nominal Height:	85mm
Diameter:	39mm
Internal Neck Diameter:	21mm
Weight:	17g
Pack Size:	72 bottles
Method of Sterilisation:	N/A (these bottles are non-sterile)
BSE/TSE Free:	Yes
Supplied with lid screwed on:	Yes
Recommended Storage Conditions:	Dry and ambient temperature

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 2-Jul-21