

1. Part Number: 8130A-XXX

2. Photograph:



3. Product Information:

Description: Translucent Square Bottles with Tamper Evident Caps

Materials of Construction: Body = High Density Polyethylene (translucent)
Cap = Polypropylene (blue)

Bottle Part Nos & Dimensions:

Part No.	Nominal Volume	Width	Capless Height	Thread Diameter	Internal Neck Diameter
8130A-100	100ml	45mm	108mm	32mm	17mm
8130A-250	250ml	60mm	136mm	32mm	17mm
8130A-500	500ml	74mm	170mm	32mm	17mm
8130A-1000	1000ml	85mm	221mm	32mm	17mm

Pressure Rating: None

Individually Bagged: No

Method of Sterilisation: N/A

Recommended Storage Conditions: Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 3-Nov-23