

1. Part Numbers: 8131A-100

2. Photograph:



3. Product Information:

Description	Sample Containers 100ml
Materials of Construction:	Body Polypropylene (PP) Lid High Density Polyethylene (HDPE)
Manufacturing Environment:	Controlled environment within factory
Nominal Volume	100ml
Brimful Volume	145ml
Colour:	Translucent with red lid
Label Area:	50 x 25mm
Graduation Marks:	Every 20ml
Overall Height (with lid on):	76mm
Internal Neck Diameter:	56mm
Internal Base Diameter:	43mm
Body Height:	72mm
Nominal Weight:	14g
Individually Bagged?	No
Pre-sterilised?	No
Number of Containers per Box:	280
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 26-Jul-24