


1. Part Numbers: 8134A-1

2. Photograph:



3. Product Information:

Description	Pasteur Pipette 1ml
Material of Construction:	Low Density Polyethylene (LDPE)
Colour:	Translucent
Graduation Marks:	Every 0.25ml
Overall Length:	155mm
Nominal Weight:	1g
Individually Bagged?	No
Method of Sterilisation:	N/A
Number of Pipettes per Box:	500
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 25-Jan-21