

1. Part Numbers: 8136H-1000A

2. Photograph:



3. Product Information:

Description	Wide Mouth HDPE Amber Bottle 1000ml	
Materials of Construction:	Bottle = HDPE	Cap = PP
	All polymers conform to:	
	• FDA 21 CFR 177.1520 • USP Class VI <87> & <88> • EU 10/2011	
Manufacturing Environment:	ISO Class 8 Cleanroom (Class 100,000)	
Nominal Volume:	1000ml	
Brimful Volume:	1011ml	
Nominal Height:	197mm	
Diameter:	91mm	
Internal Neck Diameter:	52mm	
Weight:	95g	
Pack Size:	24 bottles per pack	
Method of Sterilisation:	N/A - these bottles are non-sterile	
BSE/TSE Free:	Yes	
Supplied with lid screwed on:	Yes	
Endotoxin Free:	Yes - tested to ISO 11737-3 & USP <85>	
Particulate Free:	Yes - tested to USP <788>	
Phthalate Free:	Yes	



Leakproof:

Yes – No leaks observed when a sample bottle was filled with water, inverted with a pressure of 715mmHg for 15 minutes

Recommended Storage Conditions: Dry and ambient temperature

4. Approval:

Signed:

Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Aug-26