

**1. Part Numbers:** 8136H-500A

**2. Photograph:**



### 3. Product Information:

|                               |                                                                                                                                                   |          |
|-------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|----------|
| Description                   | Wide Mouth HDPE Amber Bottle 500ml                                                                                                                |          |
| Materials of Construction:    | Bottle = HDPE                                                                                                                                     | Cap = PP |
|                               | All polymers conform to:                                                                                                                          |          |
|                               | <ul style="list-style-type: none"> <li>• FDA 21 CFR 177.1520</li> <li>• USP Class VI &lt;87&gt; &amp; &lt;88&gt;</li> <li>• EU 10/2011</li> </ul> |          |
| Manufacturing Environment:    | ISO Class 8 Cleanroom (Class 100,000)                                                                                                             |          |
| Nominal Volume:               | 500ml                                                                                                                                             |          |
| Brimful Volume:               | 537ml                                                                                                                                             |          |
| Nominal Height:               | 165mm                                                                                                                                             |          |
| Diameter:                     | 75mm                                                                                                                                              |          |
| Internal Neck Diameter:       | 43mm                                                                                                                                              |          |
| Weight:                       | 72g                                                                                                                                               |          |
| Pack Size:                    | 48 bottles per pack                                                                                                                               |          |
| Method of Sterilisation:      | N/A (these bottles are non-sterile)                                                                                                               |          |
| BSE/TSE Free:                 | Yes                                                                                                                                               |          |
| Supplied with lid screwed on: | Yes                                                                                                                                               |          |
| Endotoxin Free:               | Yes - tested to ISO 11737-3 & USP <85>                                                                                                            |          |
| Particulate Free:             | Yes - tested to USP <788>                                                                                                                         |          |
| Phthalate Free:               | Yes                                                                                                                                               |          |



Leakproof:

Yes – No leaks observed when a sample bottle was filled with water, inverted with a pressure of 715mmHg for 15 minutes

Recommended Storage Conditions: Dry and ambient temperature

## 4. Approval:

Signed:

Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Aug-26