

1. Part Numbers: 8136H-60A

2. Photograph:



3. Product Information:

Description	Wide Mouth HDPE Amber Bottle 60ml	
Materials of Construction:	Bottle = HDPE	Cap = PP
	All polymers conform to:	
	- FDA 21 CFR 177.1520 - USP Class VI <87> & <88> - EU 10/2011	
Manufacturing Environment:	ISO Class 8 Cleanroom (Class 100,000)	
Nominal Volume:	60ml	
Brimful Volume:	67ml	
Nominal Height:	83mm	
Diameter:	39mm	
Internal Neck Diameter:	21mm	
Weight:	16g	
Pack Size:	72 bottles per pack	
Method of Sterilisation:	N/A (these bottles are non-sterile)	
BSE/TSE Free:	Yes	
Supplied with lid screwed on:	Yes	
Endotoxin Free:	Yes - tested to ISO 11737-3 & USP <85>	
Particulate Free:	Yes - tested to USP <788>	
Phthalate Free:	Yes	



Leakproof: Yes – No leaks observed when a sample bottle was filled with water, inverted with a pressure of 715mmHg for 15 minutes

Recommended Storage Conditions: Dry and ambient temperature

4. Approval:

Signed:

Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Aug-26