



1. Part Numbers: 8142P-20T

2. Photograph:



3. Product Information:

Description	Carboy 20L with Stopcock		
Material of Construction:	Body	PP	(conforms to FDA 177.1520 & USP Class VI)
	Cap	PP	(conforms to FDA 177.1520 & USP Class VI)
	Cap Seal	Silicone	(conforms to FDA 177.2600 & USP Class VI)
	Stopcock	PP	(conforms to FDA 177.1520 & USP Class VI)
Nominal Volume:	20L		
Nominal Body Height:	518mm		
Overall Height (with Cap on):	533mm		
Max. Diameter:	285mm		
Internal Neck Diameter:	65mm		
Weight:	2.0Kg		
Graduation Marks:	From 0L to 14L (every 1L)		
	Graduation marks permanently laser etched onto container body		
Autoclavable:	Yes		
Pack Size:	1		
Recommended Storage Conditions:	Dry and ambient temperature		

4. Approval:

Signed:

Name: Kuldeep Flora

Position: Quality Manager

Date: 26-Aug-26