

1. Part Numbers: 8147P-1000

2. Photograph:



3. Product Information:

Description

Media Bottle PETG - 30ml

Materials of Construction:

Bottle	Cap
- PETG - FDA 21 CFR 177.1630 - USP Class VI - BSE/TSE free - BPA free	- HDPE - FDA 21 CFR 177.1520 - USP Class VI - BSE/TSE free - BPA free

Manufacturing Environment:

ISO Class 8 Cleanroom (Class 100,000)

Nominal Volume:

1000ml

Brimful Volume:

1235ml

Nominal Bottle Height:

216mm

Width:

92 x 92mm

Internal Neck Diameter:

28mm

Weight (bottle & cap):

188g

Pack Size:

24 bottles per box (supplied as 2 packs of 12 in a single box)

Graduations:

100ml to 1000ml (every 100ml)

Method of Sterilisation:

N/A - these bottles are non-sterile



BSE/TSE Free: Yes
Supplied with lid screwed on: Yes
Endotoxin Free: Yes - tested to ISO 11737-3 & USP <85>
Particulate Free: Yes - tested to USP <788>
Leakproof: Yes - No leaks observed when a sample bottle was filled with water, inverted with a pressure of 715mmHg for 15 minutes
Recommended Storage Conditions: Dry and ambient temperature

4. Approval:

Signed:

Name: Kuldeep Flora

Position: Quality Manager

Date: 13-Oct-26