

Product Specification

Specification Number

8165A-1000

Revision

01

Part No. 8165A-1000
 Description Square Sample Bottle 1000ml
 Materials of Construction: HDPE (conforms to FDA 177.1520 & EU 10/2011)



Nominal Volume: 1000ml
 Nominal Height: 183mm
 Width: 85mm
 Internal Neck Diameter: 31mm
 Tamper-proof Closure? No
 Sterile? No
 Weight: 75g

Specification Approved By:	Signature: 	Name: Kuldeep Flora Position: Quality Manager Date: 8-Jul-19
Page 1 of 1	Sampling Systems Ltd. 4-5 Forge Mills Park, Station Road, Coleshill. Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com	