

**1. Part Number:** 8175H-150

**2. Photograph:**



**3. Product Information:**

Description:	SteriWare® Pallet Knife - 150mm
Use:	Single use, disposable Pallet Knife
Material of Construction:	High Density Polyethylene (HDPE) - virgin <ul style="list-style-type: none"> <li>• Conforms to FDA 21 CFR 177.1520</li> <li>• Conforms to EU Regulations 10/2011</li> <li>• Conforms to EC Regulations 1935/2004</li> <li>• Conforms to REACH regulation 1907/2006</li> <li>• BSE/TSE free</li> </ul>
Method of Construction:	Single piece, injection moulded
Overall Length:	260mm
Blade Length:	150mm
Blade Width:	40mm
Nominal Weight:	44g
Assembly & Primary Packing Environment:	Class 100,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	N/A
Quantity per Box:	100
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This product has not been ATEX tested

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable product. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 7-Oct-20