

**1. Part Number:** 8185P-250

**2. Photograph:**



### 3. Product Information:

Description:	GMP Tub 250ml
Use:	Wide mouth bottle with screw on cap - for use with powders
Material of Construction:	Polypropylene (PP) – virgin <ul style="list-style-type: none"> <li>• Conforms to FDA CFR 177.1520</li> <li>• Conforms to EU Regulation 10/2011</li> <li>• Conforms to EC Regulation 1935/2004</li> <li>• BSE/TSE free</li> </ul>
	White Colouring (1%)
Method of Construction:	Single piece, injection moulded
Nominal Volume:	250ml
Brimful Volume:	280ml
Method of Construction:	Single piece, injection moulded
Height:	88mm
Lid Diameter:	76mm
Nominal Weight:	56g
Supplied with Lids attached?	Yes
Individually Bagged?	No (The Tubs are packed in a box liner inside a cardboard box)
Method of Sterilisation:	N/A
Quantity per Box:	100
Recommended Storage Conditions:	Dry and ambient temperature
ATEX:	This product has not been ATEX tested

**4. Country of Origin:**

Made and packed in the UK.

**5. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 9-Apr-25