

1. Part Number: 8185P-500

2. Photograph:



3. Product Information:

Description:	GMP Tub 500ml
Use:	Wide mouth bottle with screw on cap - for use with powders
Material of Construction:	Polypropylene (PP) – virgin <ul style="list-style-type: none"> • Conforms to FDA CFR 177.1520 • Conforms to EU Regulation 10/2011 • Conforms to EC Regulation 1935/2004 • BSE/TSE free
	White Colouring (1%)
Method of Construction:	Single piece, injection moulded
Nominal Volume:	500ml
Brimful Volume:	540ml
Height:	106mm
Lid Diameter:	94mm
Nominal Weight:	83g
Supplied with Lids attached?	Yes
Individually Bagged?	No (The Tubs are packed in a box liner inside a cardboard box)
Method of Sterilisation:	N/A
Quantity per Box:	50
Recommended Storage Conditions:	Dry and ambient temperature
ATEX:	This product has not been ATEX tested

4. Country of Origin: Made in the UK.

5. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 9-Apr-25