


1. Part Number: 8220P-031

2. Label Design:



3. Product Information:

Description PharmaLabel - Sampled – Blue

Adhesive: PharmaLabel Adhesive

Colours: Blue printing on a white label

FDA Acceptable Adhesive? FDA CFR 175.105

The adhesive complies with the European food directives and legislations, FDA 175.105 and the German recommendations XXI as published by BfR. BfR (Bundesinstitut für Risikobewertung) is the German Federal Institute for Risk Assessment. The adhesive can be used in direct contact with dry and moist, non fatty foodstuffs.

The adhesive conforms to: European food regulations 1935/2004/EC, 10/2011/EU and DIN ISO 3826-1.

Dimensions: 150 x 150mm

Number of Labels per Roll: 250 (approx).

Temperature Working Temperature Range -40 to +70°C (-40 to +160°F)



Sampling Systems Ltd.

4-5 Forge Mills Park, Coleshill, Warwickshire, B46 1JH, UK.

Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com

Weight: 1300g per roll
Method of Sterilisation: N/A
Recommended Storage Conditions: Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 26-Sep-19