

**1. Part Number:** 8220X-150

**2. Photograph:**




### 3. Product Information:

Description:	Label Dispenser 150mm
Material of Construction:	<ul style="list-style-type: none"> <li>- 316L stainless steel</li> <li>- Non-slip rubber grommets fitted to base</li> <li>- HDPE rollers</li> </ul>
Max. Label Width:	Labels up to 150mm wide
Mounting:	Suitable for bench or wall mounting
Length:	350mm
Width:	162mm
Height:	100mm
Nominal Weight:	2.35kg
Recommended Storage Conditions:	Dry and ambient temperature

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 22-Apr-22