

1. Part Number: 8220X-95

2. Photograph:



3. Product Information:

Description:	Label Dispenser 95mm
Material of Construction:	- 316L Stainless Steel - Non-slip rubber grommets fitted to base - HDPE rollers
Max. Label Width:	Labels up to 95mm wide
Mounting:	Suitable for bench or wall mounting
Length:	282mm
Width:	108mm
Height:	100mm
Nominal Weight:	1.75kg
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 22-Apr-22