

1. Part Number: 8220Y-150

2. Photograph:



3. Product Information:

Description:	Label Dispenser 150mm
Material of Construction:	- Aluminium - Non-slip rubber grommets fitted to base - HDPE rollers
Max. Label Width:	Labels up to 150mm wide
Mounting:	Suitable for bench or wall mounting
Length:	350mm
Width:	162mm
Height:	100mm
Nominal Weight:	800g
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 22-Apr-22