

1. Part Numbers: 8235P-4000

2. Photograph:



3. Product Information:

Description:	Wide Mouth Square Bottle 4L
Materials of Construction:	4L Bottle Polypropylene PP (conforms to FDA 177.1520) Cap Polypropylene PP (conforms to FDA 177.1520)
Nominal Volume:	4 Litres
Nominal Height:	305mm
Width:	150mm
Internal Neck Diameter:	90mm
Weight:	530g
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 16-Aug-24