

**1. Part Number:** 8240H-1000S

**2. Photograph:**



### 3. Product Information:

Description:	SteriWare Powder Sampler - 1000mm - Sterile
Use:	Single use, disposable powder sampler
Material of Construction:	High Density Polyethylene (HDPE) - virgin <ul style="list-style-type: none"> <li>• FDA acceptable CFR 177.1520</li> <li>• EU Regulations 10/2011 compliant</li> <li>• EC Regulations 1935/2004 compliant</li> <li>• BSE/TSE free</li> </ul>
Method of Construction:	Ultrasonic welding
Nominal Length:	987mm
Diameter of Outer Tube:	40mm
Length of Sampling Slot:	60mm
Width of Sampling Slot:	10mm
Nominal Weight of Sampler:	285g
Moulding Environment:	Class 100,000 Cleanroom
Assembly & Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25kGy minimum dose)
Number of Samplers per Box:	20
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX rated

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

[change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable sampler. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 4-Apr-24