

## **Specification Sheet**

Specification Number: Revision: Date Issued: 8284R-150 04 26-Oct-20

1. Part Number:

8284R-150

2. Photograph:



## 3. Product Information:

Description:	SteriWare <sup>®</sup> Cupped Spatula Red – 150mm
Use:	Single use, powder spatula
Material of Construction:	Polystyrene – virgin, red • Conforms to FDA CFR 177.1640 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Blade Shape:	U
Blade Style:	Closed Tip
Overall Length	240mm
Blade Length:	155mm
Blade Width:	18mm
Nominal Weight:	20g
Nominal Volume (flat-filled):	21ml
Moulding & Primary Packing Environment:	Class 100,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	N/A
Number of Samplers per Box:	100
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This item has not been ATEX tested



## 4-5 Forge Mills Park, Coleshill. Warwickshire. B46 1JH. UK.

4-5 Forge Mills Park, Coleshill. Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com



**Specification Sheet Specification Number:** 8284R-150 Revision: 26-Oct-20 Date Issued:

04

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

This is a single-use, disposable sampler. Repeated use will lead to

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

7. Approval:

failure.

Signed:

Name: Kuldeep Flora Position: **Quality Manager** 26-Oct-20 Date:



Sampling Systems Ltd. 4-5 Forge Mills Park, Coleshill. Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com