

1. Part Number: 8286A-100S

2. Photograph:



3. Product Information:

Description:

SteriWare® Scoop 100ml - Sterile

Use:

Free-standing, single use, disposable scoop

Nominal Flat-filled Volume:

100ml

Material of Construction:

High Impact Polystyrene (HIPS) – rigid, white

- Conforms to FDA CFR 177.1640
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

Method of Construction:

Single piece, injection moulded

Overall Length

200mm

Handle Length:

95mm

Maximum Width:

59mm

Nominal Weight:

25g

Moulding Environment:

Class 100,000 Cleanroom

Primary Packing Environment:

Class 10,000 Cleanroom

Individually Bagged?

Yes (heat sealed PE/PET bag)

Method of Sterilisation:

Gamma Irradiation (25-45kGy)

Number of Samplers per Box:

100

Recommended Storage Conditions:

Dry and ambient temperature

Shelf Life:

5 years from the date of manufacture

ATEX:

This device has not been ATEX tested

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable scoop. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 16-Aug-24