

Specification Sheet

Specification Number: 8286A-250S Revision: 04 Date Issued: 23-Nov-22

1. Part Number: 8286A-250S

2. Photograph:





3. Product Information:

Description: SteriWare® Scoop 250ml - Sterile

Use: Free standing, single use, disposable scoop

Nominal Flat-filled Volume: 250ml

Material of Construction: High Impact Polystyrene (HIPS) – rigid, white

Conforms to FDA CFR 177.1640

• Conforms to EU Regulations 10/2011

• Conforms to EC Regulations 1935/2004

• BSE/TSE free

Method of Construction: Single piece, injection moulded

Overall Length 237mm
Handle Length: 94mm
Maximum Width: 80mm
Nominal Weight: 41g

Moulding Environment: Class 100,000 Cleanroom
Primary Packing Environment: Class 10,000 Cleanroom

Individually Bagged? Yes (heat sealed PE/PET bag)
Method of Sterilisation: Gamma Irradiation (25-45kGy)

Number of Samplers per Box: 100

Recommended Storage Conditions: Dry and ambient temperature

Shelf Life: 5 years from the date of manufacture
ATEX: This device has not been ATEX tested



Sampling Systems Ltd.



Specification Sheet

Specification Number: 8286A-250S Revision: 04 Date Issued: 23-Nov-22

4. Change Control: Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

change.control@sampling.com

5. Country of Origin: Made and packed in the UK.

6. Repeated Use: This is a single-use, disposable scoop. Repeated use will lead to

failure.

7. Approval: Signed: Name: Kuldeep Flora

Position: Quality Manager

Date: 23-Nov-22