

**1. Part Number:** 8286A-50S

**2. Photograph:**



### 3. Product Information:

**Description:**

SteriWare® Scoop 50ml - Sterile

**Use:**

Free-standing, single use, disposable scoop

**Nominal Flat-filled Volume:**

50ml

**Material of Construction:**

High Impact Polystyrene (HIPS) – rigid, white

- Conforms to FDA CFR 177.1640
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

**Method of Construction:**

Single piece, injection moulded

**Overall Length**

160mm

**Handle Length:**

75mm

**Maximum Width:**

47mm

**Nominal Weight:**

15g

**Moulding Environment:**

Class 100,000 Cleanroom

**Primary Packing Environment:**

Class 10,000 Cleanroom

**Individually Bagged?**

Yes (heat sealed PE/PET bag)

**Method of Sterilisation:**

Gamma Irradiation (25-45kGy)

**Number of Samplers per Box:**

100

**Recommended Storage Conditions:**

Dry and ambient temperature

**Shelf Life:**

5 years from the date of manufacture

**ATEX:**

This device has not been ATEX tested

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable scoop. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 23-Nov-22