

1. Part Number: 8286A-50S

2. Photograph:



3. Product Information:

| | |
|---------------------------------|--|
| Description: | SteriWare® Scoop 50ml - Sterile |
| Use: | Free-standing, single use, disposable scoop |
| Nominal Flat-filled Volume: | 50ml |
| Material of Construction: | High Impact Polystyrene (HIPS) – rigid, white <ul style="list-style-type: none"> • Conforms to FDA CFR 177.1640 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free |
| Method of Construction: | Single piece, injection moulded |
| Overall Length | 160mm |
| Handle Length: | 75mm |
| Maximum Width: | 47mm |
| Nominal Weight: | 15g |
| Moulding Environment: | Class 100,000 Cleanroom |
| Primary Packing Environment: | Class 10,000 Cleanroom |
| Individually Bagged? | Yes (heat sealed PE/PET bag) |
| Method of Sterilisation: | Gamma Irradiation (25-45kGy) |
| Number of Samplers per Box: | 100 |
| Recommended Storage Conditions: | Dry and ambient temperature |
| Shelf Life: | 5 years from the date of manufacture |
| ATEX: | This device has not been ATEX tested |

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable scoop. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 16-Aug-24