


1. Part Number: 8341H-100S

2. Photograph:



3. Product Information:

Description:	SteriWare® Long Handled Scoop 100ml - Sterile
Use:	Single use, disposable scoop
Material of Construction:	High Density Polyethylene (HDPE) - virgin <ul style="list-style-type: none"> • FDA acceptable CFR 177.1520 • EU Regulations 10/2011 compliant • EC regulation 1935/2004 compliant • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Max. Width:	50mm
Overall length:	390mm
Handle length:	300mm
Nominal Scoop Volume:	100ml
Nominal Weight:	37g
Moulding Environment:	Class 100,000 Cleanroom
Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25-45kGy)
Quantity per Box:	50
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This product has not been ATEX tested

4. Change Control: Sampling Systems operates a policy of continual improvement.



Sampling Systems Ltd.

4-5 Forge Mills Park, Coleshill, Warwickshire, B46 1JH, UK.
Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable device. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 24-Nov-22