


1. Part Number: 8352PP-100S

2. Photograph:



3. Product Information:

Description:	SteriWare® Cleaning Coupon PP - 100mm – Sterile
Use:	Single use, cleaning coupon
Material of Construction:	Polypropylene (PP) - virgin <ul style="list-style-type: none"> • Conforms to FDA 21 CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Dimensions:	100 x 100mm
Nominal Weight:	18g
Moulding Environment:	Class 100,000 Cleanroom
Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25 to 50kGy)
Quantity per Box:	50
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This product has not been ATEX tested



Sampling Systems Ltd.

4-5 Forge Mills Park, Coleshill, Warwickshire, B46 1JH, UK.
Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable product. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 24-Nov-22