


**1. Part Number:** 8352PVC-100

**2. Photograph:**



**3. Product Information:**

Description:	SteriWare® Cleaning Coupon PVC - 100mm
Use:	Single use, cleaning coupon
Material of Construction:	Polyvinyl Chloride (PVC) - virgin <ul style="list-style-type: none"> <li>• Conforms to FDA 21 Chapter 1</li> <li>• Conforms to EU Regulations 10/2011</li> </ul>
Method of Construction:	Single piece, injection moulded
Dimensions:	100 x 100mm
Nominal Weight:	28g
Moulding Environment:	Class 100,000 Cleanroom
Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	N/A
Quantity per Box:	50
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture



**Sampling Systems Ltd.**

2 Trillennium Point, Gorse Lane, Coleshill. B46 1JU. UK.

Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com

ATEX:

This product has not been ATEX tested

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

[change.control@sampling.com](mailto:change.control@sampling.com)**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable product. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 21-Aug-24