


**1. Part Number:** 8360T-50S

**2. Photograph:**



### 3. Product Information:

<b>Description:</b>	AccuPod® 1.0ml - Sterile
<b>Use:</b>	Single use, disposable sampling pod
<b>Material of Construction:</b>	High Density Polyethylene (HDPE) – virgin <ul style="list-style-type: none"> <li>• FDA acceptable CFR 177.1520</li> <li>• EU Regulations 10/2011 compliant</li> <li>• EC regulation 1935/2004 compliant</li> <li>• BSE/TSE free</li> </ul>
<b>Method of Construction:</b>	Body = Single piece, injection moulded Lid = Single piece, injection moulded
<b>Nominal Volume:</b>	5.0ml
<b>Length:</b>	61mm
<b>Width:</b>	13mm
<b>Max. Height:</b>	20mm
<b>Weight:</b>	5g
<b>Moulding Environment:</b>	Class 100,000 Cleanroom
<b>Primary Packing Environment:</b>	Class 10,000 Cleanroom
<b>Individually Bagged?</b>	Yes (heat sealed PE/PET bag)
<b>Method of Sterilisation:</b>	Gamma Irradiation (25 to 45kGy)
<b>Number of AccuPods per Box:</b>	50
<b>Recommended Storage Conditions:</b>	Dry and ambient temperature
<b>Shelf Life:</b>	5 years from the date of manufacture
<b>ATEX:</b>	This product has not been ATEX tested

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable device. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 22-Aug-24