

1. Part Number: 8370B-100S

2. Photograph:



3. Product Information:

Description:	SteriWare Funnel 100mm diameter - PP Blue - Sterile
Use:	Single use, disposable liquid funnel
Material of Construction:	Polypropylene (PP) - virgin <ul style="list-style-type: none"> • Conforms to FDA 21 CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • Conforms to REACH regulation 1907/2006 • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Max. Diameter:	117mm
ID at Top of Funnel:	100mm
Overall Height:	103mm
OD at Bottom of Spout:	11mm
ID at Bottom of Spout:	8mm
OD at Top of Spout:	17mm
Nominal Weight:	17g
Moulding Environment:	Class 100,000 Cleanroom
Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25 to 45 kGy)
Number of Samplers per Box:	50
Recommended Storage Conditions:	Dry and ambient temperature

Shelf Life: 5 years from the date of manufacture
ATEX: This product has not been ATEX tested

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable sampler. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 17-May-23