

1. Part Number: 8370H-50S

2. Photograph:



3. Product Information:

Description:

SteriWare® Funnel 50mm - Sterile

Use:

Single use, disposable funnel

Material of Construction:

High Density Polyethylene (HDPE) - virgin

- Conforms to FDA CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

Method of Construction:

Single piece, injection moulded

Max. Diameter:

68mm

ID at Top of Funnel:

50mm

Overall Height:

56mm

OD at Bottom of Spout:

8mm

ID at Bottom of Spout:

5mm

OD at Top of Spout:

12mm

Nominal Weight:

6g

Moulding & Primary Packing Environment:

Class 100,000 Cleanroom

Individually Bagged?

Yes (heat sealed PE/PET bag)

Method of Sterilisation:

Gamma Irradiation (25 to 45 kGy)

Quantity per Box:

100

Recommended Storage Conditions:

Dry and ambient temperature

Shelf Life:

5 years from the date of manufacture

ATEX:

This device has not been ATEX tested

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable device. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 21-Sep-21