

1. Part Number: 8404H-125S

2. Photograph:



3. Product Information:

Description	Wide Mouth HDPE Bottle 125ml - Sterile
Materials of Construction:	Bottle HDPE (conforms to FDA 177.1520 & USP Class VI) Cap PP (conforms to FDA 177.1520 & USP Class VI)
Manufacturing Environment:	Controlled environment within factory
Packing Environment:	Class 10,000 Cleanroom
Nominal Volume:	125ml
Brimful Volume:	140ml
Nominal Height:	99mm
Diameter:	51mm
Internal Neck Diameter:	28mm
Weight:	27g
Individually Bagged:	Yes (in a PE/PET laminate film bag)
Quantity per box:	50
Method of Sterilisation:	Gamma Irradiation (25-45 kGy)
BSE/TSE Free:	Yes
Supplied with lid screwed on:	Yes
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

Packed and sterilised in the UK.

6. Repeated Use:

This is a single-use, disposable device. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 22-Aug-24