

**1. Part Number:** 8432P-500S

**2. Photograph:**



### 3. Product Information:

**Description:**

SteriWare® Graduated Jug 500ml - Sterile

**Use:**

Single use, disposable jug

**Material of Construction:**

Polypropylene (PP) – virgin

- Conforms to FDA CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

**Method of Construction:**

Single piece, injection moulded

**Nominal Volume:**

500ml

**Brimful Volume:**

635ml

**Top Diameter (excl. spout):**

90mm

**Max Width (inc. spout and handle)**

135mm

**Base Diameter:**

75mm

**Height:**

130mm

**Nominal Weight:**

58g

**Graduation Marks:**

To 500ml (every 10ml)

**Moulding Environment:**

Class 100,000 Cleanroom

**Primary Packing Environment:**

Class 10,000 Cleanroom

**Individually Bagged?**

Yes (heat sealed PE/PET bag)

Method of Sterilisation: Gamma Irradiation (25 - 45 kGy)  
Quantity per Box: 35  
Recommended Storage Conditions: Dry and ambient temperature  
Shelf Life: 5 years from the date of manufacture  
ATEX: This device has not been ATEX tested

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

[change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable sampler. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 22-Aug-24