

1. Part Number: 8432P-500

2. Photograph:



3. Product Information:

Description:

SteriWare® Graduated Jug 500ml

Use:

Single use, disposable jug

Material of Construction:

Polypropylene (PP) – virgin

- Conforms to FDA CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

Method of Construction:

Single piece, injection moulded

Nominal Volume:

500ml

Brimful Volume:

635ml

Top Diameter (excl. spout):

90mm

Max Width (inc. spout and handle)

135mm

Base Diameter:

75mm

Height:

130mm

Nominal Weight:

58g

Graduation Marks:

To 500ml (every 10ml)

Moulding Environment:

Class 100,000 Cleanroom

Primary Packing Environment:

Class 10,000 Cleanroom

Individually Bagged?

Yes (heat sealed PE/PET bag)

Method of Sterilisation:	N/A
Quantity per Box:	35
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable sampler. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 22-Aug-24