

1. Part Number: 8555A-150S

2. Photograph:



3. Product Information:

Description:	Rounded Scissors 150mm (2 rounded blades)
Method of Construction:	Crevice free, all welds ground and polished
Material of Construction:	Stainless Steel
Surface Finish:	Satin
Overall Length:	155mm
Nominal Blade Length	55mm
Nominal Weight:	54g
Recommended Storage Conditions:	Dry and ambient temperature
Individually Bagged	Yes – Double bagged.
Method of Sterilisation	Gamma Irradiation (25 – 50 kGy)
Quantity Per Box	20
Shelf Life	2 Years from the date of packaging
Cleaning Environment	They are made in a metal working factory and then cleaned. At the packing plant they are then quality checked and wiped with IPA prior to being bagged in a class 10,000 Cleanroom.
ATEX	This device has not been ATEX tested

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 9-May-25