

**1. Part Number:** 8601S-10

**2. Photograph:**



### 3. Product Information:

Description:	Stainless Sample Splitter – 10mm slots
Material of Construction:	Sample Splitter: Stainless steel Collecting Bins: Stainless steel Frame: Painted mild steel
Method of Construction:	Spot welded
Finish:	Unpolished, welds left as laid
Number of slots:	16
Slot Width:	10mm
Overall Height:	348mm
Bin Dimensions:	273 x 137 x 84mm
Number of Bins Supplied:	2
Nominal Weight:	2.2kg (Splitter & 2 Bins)
Recommended Storage Conditions:	Dry and ambient temperature
Important Note:	Not suitable for GMP use. This product is designed for environmental use. Due to the method of construction this product is not suitable for hygienic industries (e.g. pharmaceutical and food).

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 9-Jul-24