


1. Part Number: 8601S-20

2. Photograph:



3. Product Information:

Description:	Stainless Sample Splitter – 20mm slots	
Material of Construction:	Sample Splitter:	Stainless steel
	Collecting Bins:	Stainless steel
	Frame:	Painted mild steel
Method of Construction:	Spot welded	
Finish:	Unpolished, welds left as laid	
Number of slots:	12	
Slot Width:	20mm	
Overall Height:	390mm	
Bin Dimensions:	343 x 137 x 84mm	
Number of Bins Supplied:	2	
Nominal Weight:	3.20kg (Splitter & 2 Bins)	
Recommended Storage Conditions:	Dry and ambient temperature	
Important Note:	Not suitable for GMP use. This product is designed for environmental use. Due to the method of construction this product is not suitable for hygienic industries (e.g. pharmaceutical and food).	

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 9-Jul-24