

1. Part Number: 8601S-30

2. Photograph:



3. Product Information:

Description:	Stainless Sample Splitter – 30mm slots	
Material of Construction:	Sample Splitter:	Stainless steel
	Collecting Bins:	Stainless steel
	Frame:	Painted mild steel
Method of Construction:	Spot welded	
Finish:	Unpolished, welds left as laid	
Number of slots:	12	
Slot Width:	30mm	
Overall Height:	475mm	
Bin Dimensions:	511 x 190 x 108mm	
Number of Bins Supplied:	2	
Nominal Weight:	6.6kg (Splitter & 2 Bins)	
Recommended Storage Conditions:	Dry and ambient temperature	
Important Note:	Not suitable for GMP use. This product is designed for environmental use. Due to the method of construction this product is not suitable for hygienic industries (e.g. pharmaceutical and food).	

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 9-Jul-24