

**1. Part Numbers:** 8644S6-100S

**2. Photograph:**



**3. Product Information:**

Description:	Cleaning Coupon 316Lss – 100mm – Satin Finish
Material of Construction:	316L Stainless steel
Finish:	Polished (approx. 1.2 to 0.7microns Ra)
Dimensions:	100 x 100mm
Nominal Weight:	120g
Recommended Storage Conditions:	Dry and ambient temperature

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 5-May-21