


1. Part Numbers: 8644S6-50S

2. Photograph:



3. Product Information:

Description:	Cleaning Coupon 316Lss – 50mm – Satin Finish
Material of Construction:	316L Stainless steel
Finish:	Polished (approx. 1.2 to 0.7microns Ra)
Dimensions:	50 x 50mm
Nominal Weight:	59g
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 11-May-21