

1. Part Number: 8703A-10

2. Photograph:



3. Product Information:

Description:	Titrating Automatic Burette 10ml
Product Contact Materials:	<ul style="list-style-type: none"> - Polyethylene - Borosilicate Glass - Latex
Method of Construction:	Parts push together, they can be dis-assembled for cleaning
Overall Height:	570mm
Capacity of Burette:	10ml
Graduation Marks on Burette:	Every 0.05ml (Accuracy +/- 0.03ml, Category B ISO 385)
Volume of Reservoir Bottle:	500ml
Nominal Weight:	141g
ATEX:	This device has not been ATEX tested

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 9-Jan-25