

1. Part Number: 8890P-500S

2. Photograph:



3. Product Information:

Description:	SteriWare® Powder Bottle 500ml - Sterile
Use:	Single use, disposable bottle
Material of Construction:	Polypropylene (PP) – opaque <ul style="list-style-type: none"> • Conforms to FDA CFR 177.1520 • Conforms to EU Regulation 10/2011 • Conforms to EC Regulation 1935/2004 • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Nominal Volume:	500ml
Brimful Volume:	540ml
Method of Construction:	Single piece, injection moulded
Height:	106mm
Lid Diameter:	94mm
Nominal Weight:	83g
Moulding Environment:	Class 100,000 Cleanroom
Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)

Method of Sterilisation: Gamma Irradiation (25-45kGy)
Number of Bottles per Box: 37
Recommended Storage Conditions: Dry and ambient temperature
Shelf Life: 5 years from the date of manufacture
ATEX: This product has not been ATEX tested.

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable device. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 3-Sep-24