

**1. Part Number:** 8916A-25S

**2. Photograph:**



### 3. Product Information:

<b>Description:</b>	SteriWare® Scoop with Lid 25ml - Sterile
<b>Use:</b>	Single use, disposable scoop with transparent click on lid
<b>Nominal Flat-filled Volume:</b>	25ml
<b>Material of Construction:</b>	High Impact Polystyrene (HIPS) – rigid, white <ul style="list-style-type: none"> <li>• Conforms to FDA CFR 177.1640</li> <li>• Conforms to EU Regulations 10/2011</li> <li>• Conforms to EC Regulations 1935/2004</li> <li>• BSE/TSE free</li> </ul>
<b>Method of Construction:</b>	Single piece, injection moulded
<b>Scoop Overall Length</b>	137mm
<b>Scoop Handle Length:</b>	70mm
<b>Scoop Maximum Width:</b>	38mm
<b>Nominal Weight:</b>	14g (Scoop 10g + Lid 4g)
<b>Moulding Environment:</b>	Class 100,000 Cleanroom
<b>Primary Packing Environment:</b>	Class 10,000 Cleanroom
<b>Individually Bagged?</b>	Yes (heat sealed PE/PET bag)
<b>Method of Sterilisation:</b>	Gamma Irradiation (25-45kGy)
<b>Number of Scoops per Box:</b>	50
<b>BSE/TSE Free:</b>	Yes
<b>Recommended Storage Conditions:</b>	Dry and ambient temperature
<b>Shelf Life:</b>	5 years from the date of manufacture
<b>ATEX:</b>	This device has not been ATEX tested

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable scoop. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 29-Nov-22