

**1. Part Number:** 8916A-50S

**2. Photograph:**



**3. Product Information:**

4. Description:	SteriWare® Scoop with Lid 50ml - Sterile
5. Use:	Single use, disposable scoop with transparent click on lid
Nominal Flat-filled Volume:	50ml
Material of Construction:	High Impact Polystyrene (HIPS) – rigid, white <ul style="list-style-type: none"> <li>• Conforms to FDA CFR 177.1640</li> <li>• Conforms to EU Regulations 10/2011</li> <li>• Conforms to EC Regulations 1935/2004</li> <li>• BSE/TSE free</li> </ul>
Method of Construction:	Single piece, injection moulded
Scoop Overall Length	160mm
Scoop Handle Length:	75mm
Scoop Maximum Width:	47mm
Nominal Weight:	23g (Scoop 14g + Lid 9g)
Moulding Environment:	Class 100,000 Cleanroom
Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25-45kGy)
Number of Samplers per Box:	50
BSE/TSE Free:	Yes
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested

**6. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**7. Country of Origin:**

Made and packed in the UK.

**8. Repeated Use:**

This is a single-use, disposable scoop. Repeated use will lead to failure.

**9. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 29-Nov-22